

PMT Non-Linearity Pre-Tests Checklist

PMT Serial:..... Name:.....(Initials)

Start date:..... (YYYY-MM-DD) Time:.....(hh-mm)

Base-stages:..... High Voltage:.....V Operational Time:..... h

Notes:.....

01. Stop the data collection script ----- ☐
02. Turn off all the power supplies (**HV, DCsupply, Photodiode Temp., Photodiode Amp.**) ☐
03. Note down the operational time for the previous PMT (**on paper and spreadsheet**) ----- ☐
04. Swap the PMT ----- ☐
05. Connect the **unitary gain base** to the new PMT **with HV cable** ----- ☐
06. Remove the PMT cap and set the distance between PMT and filter wheel to **14cm** ----- ☐
07. Cover the PMT and filter wheel with the cardboard dark box ----- ☐
08. Lock the metal dark box and place the key on the designated location ----- ☐
09. Set HV=**280V**, DC-ch3=**2V (Constant LED)**, Filter-Position=**9** ----- ☐
10. Record the currents using the Picoammeter (**Wait until you read 0.01nA on the meter**) ☐

For 4 stage			For 3 stage		
	Constant LED Voltage (V)	Cathode Current (nA)		Constant LED Voltage (V)	Cathode Current (nA)
-v7			7 nA	-v12	12 nA
-v9			9 nA	-v15	15 nA
-v12			12 nA	-v18	18 nA

11. Unplug the Picoammeter from the power outlet (to reduce AC noise) ----- ☐
12. Turn off the HV power supply ----- ☐
13. Connect **3/4 stage base** to the new PMT **with HV cable** ----- ☐
14. Adjust the distance between PMT and Filter wheel to be **14cm** ----- ☐
15. Cover the PMT and filter wheel with the cardboard dark box ----- ☐
16. Lock the dark box and place the key on the designated location ----- ☐
17. Turn on High Voltage Power Supply (**4-Stage base: -600V, 3-stage base: -800V**) ----- ☐
18. Turn on 3 Channel DC Power supply **[I/O]** ----- ☐
19. Turn on the photodiode temperature controller ----- ☐
20. Plug in the photodiode amplifier barrel connector ----- ☐
21. Run the “**record.sh**” script ----- ☐
22. Wait for the first run to complete ----- ☐
23. Adjust the High Voltage if the max anode current is not within the **8-10uA** range. Input the new high voltage into the data collection script. (**HV_{min}= -500V, Hv_{max}= -900V**) ----- ☐

Note: If anode current cannot be set to the required range by adjusting HV, change the base